

OPERATING THE CF-1000BRL HIGH VOLUME SAMPLER

Purpose This Air Quality Group procedure describes the steps to operate the CF-1000BRL High Flow sampler used for sampling air for the AIRNET Project.

Scope This procedure applies to the individuals assigned to operate and maintain the air sampler.

In this procedure This procedure addresses the following major topics:

Topic	See Page
General Information About This Procedure	2
Who Requires Training to This Procedure?	2
Worker Safety	4
Setting Up The Sampler	5
Calibration	7
Records Resulting from This Procedure	8

Hazard Control Plan The hazard evaluation associated with this work is documented in Attachment 1: Initial risk = **low**. Residual risk = **low**. Work permits required: none. First authorization review date is one year from group leader signature below; subsequent authorizations are on file in group office.

Signatures

Prepared by: Daniel Herrera, ESH-17	Date: <u>02/13/01</u>
Approved by: Craig Eberhart, Air Quality Monitoring Project Leader	Date: <u>2/13/2001</u>
Approved by: Terry Morgan, QA Officer	Date: <u>2/13/2001</u>
Work authorized by: Doug Stavert, ESH-17 Group Leader	Date: <u>2/14/2001</u>

3/6/2001

CONTROLLED DOCUMENT

This copy is uncontrolled if no signatures are present or if the copy number stamp is black. Users are responsible for ensuring they work to the latest approved revision.

General information about this procedure

Attachments This procedure has the following attachments:

Number	Attachment Title	No. of pages
1	Hazard Control Plan	2

History of revision This table lists the revision history and effective dates of this procedure.

Revision	Date	Description Of Changes
0	2/20/01	New document.

Who requires training to this procedure? The following personnel require training before implementing this procedure:

- Personnel assigned to operate and change samples on the CF-1000BRL sampler.

Training method The training method for this procedure is **on-the-job** training by a previously-trained individual and is documented in accordance with the procedure for training (ESH-17-024).

Prerequisites In addition to training to this procedure, the following training is also required prior to performing this procedure:

- First Aid
- Cardiopulmonary Resuscitation (CPR)
- Field Survival
- ESH-17-202, "Environmental Sampling of Ambient Airborne Radionuclides"
- ESH-17-207, "Evaluation of AIRNET Sampler Sites Against Siting Criteria" (if determining sampler location)

General information, continued

**Definitions
specific to this
procedure**

None.

References

The following documents are referenced in this procedure:

- ESH-17-024, “Personnel Training”
 - ESH-17-202, “Environmental Sampling of Ambient Airborne Radionuclides”
 - ESH-17-207, “Evaluation of AIRNET Sampler Sites Against Siting Criteria”
 - Operating Manual for CF-1000BRL High Volume Sampler
-

Note

Actions specified within this procedure, unless preceded with “should” or “may,” are to be considered mandatory guidance (i.e., “shall”).

Worker Safety

Performing work safely

DO NOT perform work under conditions you consider unsafe. Before beginning work described in this procedure, review safety needs and requirements, identify hazards, and develop hazard mitigation measures. Be aware that facility configurations and hazards may change between visits. Hazards to assess include, but are not limited to the following:

Facility management units - Work control is the responsibility of the Facility Manager in whose area one may want to locate a sampler. Obtain approval from facility management before beginning work to locate a monitor in a Facility Management Unit. Facility management must also have knowledge of your presence and activity during subsequent normal operations. Ensure you have completed all facility-specific training requirements.

Contact your supervisor and the project leader if working conditions are found to be unsafe.

Setting up the sampler

Description of system The CF-1000BRL is a high flow sampler is Portable, Maintenance-free, High Volume Air Sampling System that can be used for continuous or intermittent sampling. The CF-1000 is a 110VAC powered sampler that will pull air at any rate up to 50 CFM. The particles can be collected on most any desired filter media. The AIRNET project uses polypropylene filter media (same as used to collect AIRNET samples according to ESH-17-202) and samples at a rate of 40 cfm. The CF-1000 comes standard with a 4" diameter filter holder. The sampler can be adapted from the 4" diameter ring to an 8" x 10" holder by simply threading on a (HI-Q model#) CFPH-810.

Siting the sampler The sampler will generally be used in remote areas where a generator will provide power if utility power is not available. The **Air Quality Monitoring Project Leader** will determine the location.

Guidance for siting the station is provided in ESH-17-207 ("Evaluation of AIRNET Sampler Sites Against Siting Criteria").

Use a GPS unit to collect site location coordinates.

Enter special studies info and get c-of-c form Ensure the data about the sampling event is entered into the AIRNET Access database under Special Studies. Obtain the chain-of-custody form for the sampling event from this database. Contact the database owner (analytical chemistry coordinator) for assistance.

Set up the sampler Set up the sampler unit so that the air inlet is free from obstructions. For most sampling, the sampler will run in a constant flow, where the sampler will run continuously and the samples will be collected as determined by the Air Monitoring Project Leader.

Steps to set up the CF-1000BRL To set up the CF-1000BRL for a sampling job, perform the following steps:

Step	Action
1	Set up tripod, extend legs to maximum height. Be sure to set firmly in the ground or use sand bags to secure it.

Steps continued on next page.

Setting up the sampler, continued

Step	Action
2	Remove the wing nut from fixed bolt located on the bottom of the CF-1000BRL. Place the sampler on the tripod and replace wing nut, tighten firmly.
3	Attach the filter holder by threading it on to the sampler.
4	Loosen the four corner knurled nuts from the filter holder. Remove the top half of the filter holder.
5	Place filter on to holder and reinstall the top half of the filter holder. Tighten the four knurled nuts firmly (do not over tighten.)
6	Connect sampler to GFCI-equipped power source or use a GFCI extension cord. Double check connections before energizing equipment. Use only "W"(wet) rated extension cords.
7	Turn the power switch on.
8	Adjust the flow rate control knob to 40 cfm.
9	Reset the sample timer to begin sample.
10	Record necessary start information on the chain-of-custody form.
11	When sample time is complete, turn off the sampler and unplug the power cord.
12	Complete the chain-of-custody form to document the sample collection.
13	Remove filter and place it in its properly identified poly bag.

Submitting sample for analysis

Follow the instructions of the Air Monitoring Project Leader for submitting the collected filter samples. Include the c-of-c form with the sample and keep a copy in the "AIRNET Special Studies Field Data Validation And Verification Notebook".

Maintenance

The CF-1000BRL High Flow Sampler is maintenance-free.

Calibration

Re-calibrating the sampler before each use Before each use, calibrate the sampler.

Steps to calibrate the CF-1000BRL

To calibrate the sampler, perform the following steps:

Step	Action
1	Install filter media, new filter holder, and filter holder adapter into the pump intake. IMPORTANT: Ensure you are using the same filter media as will be used for sample collection.
2	Connect Hi-Q air flow calibrator to the filter holder adapter.
3	Unplug power cord from outlet.
4	Remove the front panel screws and separate the front panel from the housing to expose the rotameter needle valve which has a red dust cover over it. Remove the dust cover.
5	Support the face panel securely in a vertical position for proper operation of the rotameter.
6	Plug in power and turn on power.
7	Adjust flow rate for a calibrator flow rate near the normal operating rate.
8	Read the flow rate on the rotameter. It should be the same as the calibrator flow rate within 5% full scale. If not: Re-calibrate flow meter by adjusting the small brass needle valve located near the top of the rotameter. Adjust until the rotameter indicates the same flow rate as the calibrator. If matching the flow rate is not possible it may be necessary to re-mark the existing scale or purchase a blank scale and re-mark it.

Records resulting from this procedure

Records

The following records are generated as a result of this procedure:

- entries in HiVol Sampler logbook (logbook will be submitted according to ESH-17-011)
- chain-of-custody form for sample (shipped with sample and returned with sample results; copy in “AIRNET Special Studies Field Data Validation And Verification Notebook”)

HAZARD CONTROL PLAN

1. The work to be performed is described in this procedure.

“Operating the CF-1000BRL High Volume Sampler”

2. Describe potential hazards associated with the work (use continuation page if needed).

Lifting Injuries--lifting pumps in/out of air sampling station; lifting/moving new housing and materials (timbers).

Abrasions and bruises from moving pumps, housing, using hand tools, and moving/handling materials (timbers)

Falls/tripping

Animal Injuries (snakes, spiders, mountain lions, etc.)

Weather (lightning)

High Explosives testing (TA-15, TA-16, TA-49, etc.)

Radiation Areas (TA-54- Area-G, TA-16, etc.)

Electrical shock in wet conditions

Dropping materials on feet

3. For each hazard, list the likelihood and severity, and the resulting initial risk level (before any work controls are applied, as determined according to LIR300-00-01.0, section 7.2)

Abrasions /bruises from moving pumps, housings and materials--Occasional/Negligible = Minimal

Lifting Injuries--lifting pumps and housings --Occasional/Moderate = Low

Falls/Tripping - Moderate/Occasional = Low

Animal Injuries - Critical/Remote = Minimal

Lightning - Catastrophic/Remote = Low

Entry into High Explosives testing Areas -- Critical/Remote = Minimal

Entry into posted Radiation/Controlled areas -- Negligible/Remote = Minimal.

Electrical shock in wet conditions -- Catastrophic/Remote = Low

Dropping materials onto feet -- Critical/Improbable = Low

Overall *initial* risk: ☐ Minimal ☒ Low ☐ Medium ☐ High

4. Applicable Laboratory, facility, or activity operational requirements directly related to the work:

☐ None ☒ List:

Work Permits required? ☒ No ☐ List:

LIR-402-706-01 “Personnel Dosimetry”

LIR-402-718-01 “Radiological Training”

Access Control Requirements for TA-15, TA-16, TA-49, TA-54

29CFR1926.500, Subpart M, Section 502, “Fall protection”

National Fire Protection Code—for use of electrical GCFIs

LIR 402-600-01.0 “Electrical Safety” for all electrical hazards

HAZARD CONTROL PLAN, continued

5. Describe how the hazards listed above will be mitigated (e.g., safety equipment, administrative controls, etc.):

Abrasions/bruises from moving pumps and materials -- Use common sense to avoid these injuries; also covered under Employee Orientation training

Lifting Injuries -- lifting pumps and housings -- training in proper lifting techniques.

Falls/tripping -- the ESH-17 Employee Orientation includes training and awareness of tripping and falls.

Animal Injuries -- same as above.

Lightning -- same as above.

Entry into High Explosives testing areas -- existing controls are stringent and not easily bypassed.

Entry into posted Radiation/Controlled areas -- Area-G and TA-15 controls are stringent and not easily bypassed.

Electrical shock in wet conditions -- all stations have been retrofitted with GFCI (ground fault circuit interrupts)

Dropping materials onto feet -- Steel-toed shoes or boots are required anytime pumps, station houses, timbers, or other heavy equipment is moved.

6. Knowledge, skills, abilities, and training necessary to safely perform this work (check one or both):



Group-level orientation (per ESH-17-032) and training to this procedure.



Other → See training prerequisites on procedure page 3. Any additional describe here:

7. Any wastes and/or residual materials? (check one) ☒ None ☐ List:

(Used filter media or analytical residue is disposed by the analytical laboratory.)

8. Considering the administrative and engineering controls to be used, the *residual* risk level (as determined according to LIR300-00-01.0, section 7.3.3) is (check one):



Minimal



Low



Medium (requires approval by Division Director)

9. Emergency actions to take in event of control failures or abnormal operation (check one):



None



List:

For all lifting, trips, falls, burns, cuts, electrical shocks, and animal related injuries, provide first aid and see that injured person is taken to ESH-2 or the hospital. For any exposed, energized electrical wires, contact JCI or the appropriate authority to turn off the power. Follow all site specific emergency plans for any radiation or explosives emergencies.

Signature of preparer of this HCP: This HCP was prepared by a knowledgeable individual and reviewed in accordance with requirements in LIR 300-00-01 and LIR 300-00-02.

Preparer(s) signature(s)

Name(s) (print)

/Position

Date

Signature by group leader on procedure title page signifies authorization to perform work for personnel properly trained to this procedure. This authorization will be renewed annually and documented in ESH-17 records.

Controlled copies are considered authorized. Work will be performed to controlled copies only. This plan and procedure will be revised according to ESH-17-022 and distributed according to ESH-17-030.